



## How to visualize (the dynamics of) molecular interactions?

### Förster or Fluorescence Resonance Energy Transfer (FRET)

#### Intramolecular FRET (biosensors)

##### METHOD

Ratiometric imaging approach  
(Donor/Acceptor)

##### MICROSCOPES

- confocal: Nipkow/ spinning disk, Nikon A1r spectral/CLSM
- widefield: Bioflux

##### SOFTWARES

NIS-Elements, Metafluor, Metamorph, ImageJ

##### CONDITIONS

- ❗ Key factors are:
  - Distance A-D  $< 100 \text{ \AA}$  ( $R_0$ )
  - Overlap  $\lambda_{em}(\text{Donor}) - \lambda_{ex}(\text{Acceptor})$
  - Fluorophore orientation
    - Biological context
    - Transfection level
- ❗ Controls are required to define the FRET efficiency (lowest and highest values, FRET range):
  - positive control (high FRET efficiency, biosensor constitutively active)
  - negative control (low FRET efficiency, biosensor constitutively inactive)

##### PROTOCOLE

Kardash et al. (2011). Imaging protein activity in live embryo using fluorescence resonance energy transfer biosensors. Nat. Protoc. Nov 3;6(12):1835-46

<http://zeiss-campus.magnet.fsu.edu/articles/spectralimaging/spectralfret.html>

#### Intermolecular FRET

##### METHODS

- FRET by intensity measurement
  - FRET by Acceptor photobleaching
  - FRET by Donor photobleaching
  - FRET by lifetime measurement (TCSPC, Time Correlated Single Photon Counting)
  - FRET by fluorescence anisotropy

##### MICROSCOPES

- Confocal: Nikon A1r spectral/CLSM
- FLIM (Fluorescence Lifetime Imaging Microscopy) by single-photon counting
- FCS/ FCCS (Fluorescence (Cross-) Correlation Spectroscopy)

##### SOFTWARES

NIS-Elements, ImageJ

##### CONDITIONS

- ❗ Key factors are:
  - Distance A-D  $< 100 \text{ \AA}$  ( $R_0$ )
  - Overlap  $\lambda_{em}(\text{Donor}) - \lambda_{ex}(\text{Acceptor})$
  - Stoichiometry D-A 1:1
  - Fluorophore orientation
  - Biological context
  - Transfection level
- ❗ Controls are required to define the FRET efficiency (lowest and highest values):
  - D alone, A alone, D+A (constructs)

##### PROTOCOLE

<https://www.unige.ch/medecine/bioimaging/en/information/tutorials/> and go to: "F-techniques" \* How to perform FRET experiments?"

Broussard et al. (2013). Fluorescence resonance energy transfer microscopy as demonstrated by measuring the activation of the serine/threonine kinase Akt. Nat. Protoc. 8, Jan 10, 265-281.

