

TIRF


Total Internal Reflection Fluorescence microscopy or Evanescent Wave Microscopy

■ Live imaging and observations of fixed samples

- Green and Red fluorescence

■ Temperature, Humidity, CO2 control

■ Observations of biological structures/ proteins near the plasma membrane



Some key TIRF formulas

! **Snell's law:**
 $n_1 \sin \theta_1 = n_2 \sin \theta_2$
 (θ : angle of incidence, n : refractive index)



! **Critical angle:**
 at the critical angle θ_c , $\theta_2 = 90^\circ$; $\sin 90^\circ = 1$;
 $n_1 \sin \theta_c = n_2$
 $\theta_c = \sin^{-1}(n_2/n_1)$
 If $n_1 = 1.515$ and $n_2 = 1.36$, $\theta_c = \sin^{-1}(1.36/1.515) = 63.85^\circ$

! **Evanescent field:**


$$I_z = I_0 \exp^{-z/d_p}$$

$$d_p = \frac{\lambda}{4\pi \sqrt{n_1^2 \sin^2 \theta_1 - n_2^2}}$$

I : intensity
 z : distance
 λ : wavelength
 d_p : penetration of depth



Bleaching	★★★★☆	low
Resolution	★★★★☆	medium
Speed acq.	★★★★☆	good
Multicolor	★★★★☆	3
Depth (z)	★★★★☆	<100 nm



Lasers 1

405 nm 50 mW
 488 nm 50 mW
 561 nm 25 mW

Mounted filters 2

TIRF Quad Filterset 405/488/561/640
 TIRF FITC : ex 488/10 DM502 em 530/43
 TIRF TRITC : ex 561/10 DM580 em 605/50

Objectives 3

10x 0.3 Plan Fluor DIC C / N1 WD 16mm
 40x 1.3 Oil Plan Fluor DIC H WD 0.2mm
 63x 1.4 Oil Plan Apo DIC H WD 0.21mm
 100x 1.49 Oil CFI Apochromat TIRF

Detector 4

Camera Andor iXon+
 (Model : DU-897E-CS0-#BV-500)
 Pixels : 512 x 512 (16um x 16um)
 56 fps full frame
 Cooling -80°C

Microscope 5

Nikon Eclipse Ti with PFS

Software 6

- Windows 7
 - NIS elements