

Recommendations for STED nanoscopy: 5 key points for sample preparation

- 1 As secondary antibodies (optimized for STED imaging), please use one or several shown in the list below:
[Antibodies surrounded in red are available at the Bioimaging, contact us for more information]

Secondary Antibodies	Alternative to	Excitation	Depletion
Alexa Fluor 488	ATTO488 or Dylight 488	@488 nm	@ 592 nm
Abberior STAR 488	ATTO488, AF488, Dylight 488	@488 nm	@ 592 nm
AlexaFluor 568, 555, 532		from 532 to 568 nm	@ 592 nm
Abberior STAR 580	ATTO590/594, AF 584/594	from 550 to 590 nm	@ 775 nm
Alexa Fluor 594	ATTO590/594, STAR 580	from 550 to 590 nm	@ 775 nm
Alexa Fluor 647	STAR RED, ATTO647N, STAR 635P, AF647, Cy5	diode lasers @635 or 650 or @647 with a Krypton laser	@ 775 nm
Abberior STAR 635P	STAR RED, ATTO647N, AF647, Cy5	diode lasers @635 or 650 or @647 with a Krypton laser	@ 775 nm
Abberior STAR RED	ATTO647N, AF647 or Cy5	diode lasers @635 or 650 or @647 with a Krypton laser	@ 775 nm

*** Examples of Dyes combinations: STAR580 +STAR635P (2 colors), AF488 + AF594 (2 colors), STAR580 + STAR RED + AF488 (3 colors), AF488 + AF594 + STAR RED (3 colors), etc....***

*** Good results have been reported with the following fluorescent proteins: mTurquoise2, mTFP1, eGFP, EmGFP, eYFP, Venus, mCitrine, DsRed/ mRFP, mStrawberry, mCherry...***




- 2 Do not use “Blue dyes” such as DAPI, Hoechst or even Propidium ionide or ethidium Bromide for labeling or in embedding media. These dyes may be excited by the STED laser and thereby might lead to high background and blurred images. As a good alternative, please use PICOGREEN (available at the Bioimaging), YOYO-3 or TO-PRO-3 (InVitrogen) to label the nucleus.
- 3 Use one of the following mounting medium: Glycerol/ PBS with a fresh antifading (NPG or Dabco), Prolong Antifade Gold (available at the Bioimaging) or Diamond, Fluoromount G, Mowiol, Abberior Liquid Antifade and Solid Antifade, Roti©-Mount FluorCare.

*** Please do not use Vectashield, Vectashield Hard set or other embedding media containing p-phenylenediamine as antifading reagent.***

- 4 Choose the correct coverglass: use the 1.5 coverslips (0.170 +/- 0.01 mm thick, Hecht-Assistent, cat. Number: 1014/ 2424) (available at the Bioimaging).
- 5 Seal the coverslips with nail polish.

The Bioimaging core facility can provide you a kit - if necessary- which contains the following components:

- secondary Ab (anti-rat, anti-mouse or anti-rabbit): Abberior STAR RED and Abberior STAR 580
- PicoGreen stock solution
- Coverslips
- Aliquots of mounting medium (Prolong Antifade Gold)



More information about STED principle at :

* <https://www.leica-microsystems.com/science-lab/sted-super-resolution-microscopy-nanoscopy-principles-and-photophysics/>

* <http://www.abberior.com/knowledge/microscopy-tutorials/sted/>

Bioimaging
CORE FACILITY

Faculty of medicine
Geneva