

Dear NMR users,

A new composite experiment is now available on the 400 MHz named **MP\_NOAH3\_SBC**.

It is similar to the “**FULL\_ANALYSIS**” composite experiment but in this case, 2D experiments (HSQCed, HMBC and COSY) are combined into an NMR “supersequence” (**NOAH – NMR by Ordered Acquisition using 1H-detection**)<sup>1</sup> that can be recorded in a single measurement involving a single recovery delay.

The data collection time is considerably reduced (25 minutes). Usually when running the individual experiment COSY, HSQC and HMBC experiments with 2 scans, the total experimental time is of 45 minutes.

Below are instructions for use of the supersequence NOAH in iconNMR:

A new composite experiment was implemented on both 400 and 500 MHz spectrometers:  
**MP\_FULL\_NOAH3\_SBC**


Once having selected the **MP\_NOAH3\_SBC** composite within IconNMR, the following three experiments (1H, 13C and NOAH) will appear:

Holder	Type	Status	Disk	Name	No.	Solvent	Experiment
▼ 1	3	Available					
		Available	/opt/nmrdata/service	MP_NOAH	10	CDCl3 chloroform-d	N PROTON D 1D 1H Normal (NS=16)
		Available	/opt/nmrdata/service	MP_NOAH	11	CDCl3 chloroform-d	N MP_zgpg30
		Available	/opt/nmrdata/service	MP_NOAH	12	CDCl3 chloroform-d	N MP_NOAH3_SBC
	F2		/opt/nmrdata/service	MP_NOAH	10		

For **MP\_NOAH3\_SBC** HSQC, HMBC and COSY experiments will be run **at once within 25 mins**. At the end of the acquisition, an automatic script splits the single NOAH data set into 3 individually processed 2D data sets.

For example, if the NOAH experiment is acquired in EXPNO 12, it will produce data sets numbered:

- 12001 (HSQC)
- 12002 (HMBC)
- 12003 (COSY)



<sup>1</sup> E. Kupce, T. D. W. Claridge, Chem. Commun., 2018, 54, 7139.